

# DETECTION TIMES

**APPROXIMATE\*** AMOUNT OF TIME DETECTABLE IN ORAL FLUID

DRUG	ABBREVIATION	DETECTION PERIOD
AMPHETAMINE	AMP	24 - 36 HOURS
BENZODIAZEPINE	BZO	24 - 36 HOURS
COCAINE	COC	24 - 36 HOURS
CANNABIS	THC	< 24 HOURS
METHADONE	MTD	24 - 36 HOURS
METHAMPHETAMINE	MET	24 - 36 HOURS
OPIATES	OPI	24 - 36 HOURS
PHENCYCLIDINE	PCP	24 - 36 HOURS

For further help or information call Quantum Diagnostics on 01992 651 111

*\* TO BE USED AS AN APPROXIMATE GUIDE ONLY.*